

HYBRIDSELECT®

Mattress Overlay



User-Friendly Control Unit

Rapid inflation, easy set-up and a motor so quiet, the patient will hardly know it's there!

The SelectAir® System Quietly circulates air through individual cells and allows air to escape from the mattress.

Wicking Action: Air escapes through tiny holes on the top of the mattress and circulates under a vapor-permeable cover to wick moisture away from patient.

HYBRIDSELECT® Features

HealFloat: Supple-Soft Foam heel cushion provides added protection in the foot area of the mattress and helps to reduce pressure points on the heels.

Pulsation (alternating pressure): Provides two cycle settings that continually readjust the system to ensure patient comfort.

Quick Setup: HybridSelect sets up in minutes on top of the existing mattress. No need to store the mattress!

Mattress overlay with zoned, "true" low air loss used in acute care hospitals nationwide with Supple-Soft Foam heel cushion in foot section.

Weight Limit: 300 lbs. (136 kg)

Mattress Overlay: 80" x 36" x 5.5"

Cover: Water resistant, vapor-permeable, flame resistant and biocompatible

Warranty: One Year Limited Warranty on control unit and soft goods;

Medicare Code: E0372 Group 2 Support Surface

LIT-0004 | R010422


Manufactured in the USA by
MOXI
ENTERPRISES
Innovative health solutions

1939 S. Vandeventer Ave.
St. Louis, MO 63110

T 1-800-942-7678
E info@moxiusa.com
W moxiusa.com