

SELECTPROTECT

SAFETY ENHANCEMENT SYSTEM



Concerned about patient falls and bed safety? Afraid of questionable aftermarket add-ons?

Provide real patient safety and protection. Introducing **SelectProtect**[®], the simple safety enhancement system that helps reduce the risk of patient falls and promotes bed-bound safety.

SelectProtect[®] is designed to be used exclusively with standard MOXI **SelectAir**[®] or **SelectAir Max**[®] **Low Air Loss Mattress Replacement Systems**, and without bed side

rails*. The unique built-in bolsters help to promote a safe sleeping environment with compromising the patient's sleeping area.

Don't risk patient safety with aftermarket add-ons that are not designed to fit standard **SelectAir**[®] low air loss mattress systems. Choose **SelectProtect**[®].

* Use of **SelectProtect**[®] should be determined based upon facility protocols, applicable laws, regulations, guidelines and individual user safety considerations. It does not replace the need for patient fall prevention measures.



REAL PROTECTION & SAFETY

Patient safety is always a major concern. *SelectProtect*[®] is designed to be used exclusively with standard MOXI *SelectAir*[®] or *SelectAir*[®] Max Low Air Loss Mattress Replacement Systems to increase patient bed safety. Integrated side bolsters help provide a safe sleeping environment for patients.

SelectProtect[®] is designed to be used with bed side rails¹. It is ideal for use in many healthcare settings where bed side rails are not used including long-term care facilities and in the home.

COMPATIBILITY

SelectProtect[®] is specifically designed to fit standard *SelectAir*[®] and *SelectAir*[®] Max Low Air Loss Systems. It replaces the standard cover and easily zips around the *SelectAir*[®] mattress. Simply order *SelectProtect*[®] as an option with new systems or as an upgrade to systems that are already in use.

REDUCE RISK

Independent testing has demonstrated that *SelectProtect*'s foam bolsters help prevent subjects from falling out of a low air loss mattress replacement system². Choose *SelectProtect*[®] to promote bed-bound safety.

APPLICATIONS

The *SelectProtect*[®] system is most effective for users who:

- are at risk of falling out of bed
- reside in facilities where bed side rails* are not used
- currently use standard *SelectAir*[®] or *SelectAir Max*[®] low air loss mattress systems and need a safe sleeping environment

FEATURES & BENEFITS

Protection - Integrated side bolster system helps protect bed-bound patients from falls.

Safety - Built-in bolsters help promote a safe sleeping environment for the patient while shoulder bolsters do not compromise sleeping area.

Compatible - Designed to fit standard *SelectAir*[®] or *SelectAir Max*[®] low air loss mattress systems already in use.

Durable - *SelectProtect*'s vapor-permeable, fluid-resistant material is durable and easy to clean.

SPECIFICATIONS

Includes - *SelectProtect*[®] cover, 2 shoulder bolsters, 2 foot bolsters, operation instructions.

Construction:

Top: 100% polyurethane coated polyester

Bottom: PVC-backed polyester (DEHP and BPA free).

Flame resistant: Meets CFR 16 Part 1632.

Size - For use with standard size *SelectAir*[®] or *SelectAir Max*[®] low air loss mattress systems only (36 inch mattress size). Custom sizes are not available.

Weight Limit - *SelectProtect*[®] is intended for patients that weigh 75 to 400lbs. which is the same weight limit for standard *SelectAir* and *SelectAir* low air loss mattress replacement systems.

Warranty - 1 Year Limited Warranty

¹ Use of *SelectProtect* should be determined based upon facility protocols, applicable laws, regulations, guidelines and individual user safety considerations. It does not replace the need for patient fall prevention measures.

² Independent testing was conducted using *SelectProtect* as a cover replacement for a standard *SelectAir* low air loss mattress replacement system and using a mannequin as the patient.